

Tek – Air FVC 2000 Plus

FVC 2000 Plus (Face Velocity Control)

This fume hood controller display consists of the liquid crystal display, LED indicators and three momentary push button switches.

Mute – silences alarm and LED remains on until condition changes.

Pressing and holding the Mute button when there is no alarm will cause the display to perform a self test.

Emergency Purge - maximum opening of hood exhaust valve activated - LED red on

Alarm – red LED flashes in emergency until mute pressed (as above)

Normal – green LED on

Parameters – displays programmed air flow values and set points (used in conjunction with FVC 2600 Fume hood con-



Closing the Sash can save \$1500 and 50,000 lbs of CO₂ per year.

QUESTIONS?

If your chemical fume hood needs servicing due to low or no air flow, broken dome lights, broken sash (or sash is difficult to slide up or down), air flow sensor malfunction, or the hood needs annual re-certification, please fill out the chemical fume hood service request at (<https://www.case.edu/finadmin/does/LabSafety/fume.html>) or contact

Mary Ellen Scott, PhD
216.368.6077

ENVIRONMENTAL HEALTH AND SAFETY

Case Western Reserve University
Service Building 1st Floor
2220 Circle Drive, First Floor
Cleveland, Ohio 44106-7227

Office: 216.369.2907
Fax: 216.368.2236
Web: [does.case.edu](https://www.case.edu)

ENVIRONMENTAL HEALTH AND SAFETY

FUME HOOD MONITORS II



Monitor Mania

We use many different kinds of air flow monitors for our laboratory fume hoods. All of these however have a singular purpose: TO INDICATE AIR FLOW INSIDE THE HOOD. These monitors are calibrated and an alarm will sound if the air flow is above or below the a safe set point as determined by experiment, hood type and room air.



TeK- Air: Apex Premier

Green light on – normal, flashes for warning condition



Red light on – indicates alarm condition, flashes on mute

Yellow light on- indicates purge button has been depressed

Mute button – silences alarm, programmable active/inactive mode for changing flow set points. Re-beep mode – reminder that alarm is muted.

Purge button – activates purge relay, programmable

Scroll – enables user to review entire monitor configuration (Basic models do not have digital readout)



TEL: AFA 1000

The operator display panel features:

Velocity bar graph or **Alarm Time Line**

Velocity digital display

LED Indicators (green - normal, yellow-caution and red - low)

Enter button - momentarily mutes alarm and after pressing (5 sec) gives access to calibration and configuration menus using (+ / -) scrolling. (Basic models have an air flow sensor, LED indicators.

Set button - calibration. (The Kewanee Air Alert 600 is identical to the AFA 1000. Some



TEL: AFA 500



Flow Safe: VFV Fume Hood Alarm and Baffle Control

VFV (variable face velocity) controls: commonly used on Vortex II fume hood.

The VFV fume hood alarm monitors fume hood performance by monitoring the strength of the bi-stable vortex.

Led Indicators: green (normal), yellow (caution) red (alarm—low flow)

Mute Button - silences alarm and the **Audio off** light illuminates. Alarm rearms once condition clears.

To test monitor, cover the sensing hole (port with plastic nut) located on inside of wall of hood, alarm sounds. (Some Air Sentry Fume



Air Sentry